

Product datasheet

Anti-Striatin antibody ab123291

1 图像

概述

产品名称	Anti-Striatin抗体
描述	兔多克隆抗体to Striatin
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 172-202 of Human Striatin (NP_003153.2).
阳性对照	NCI-H292 cell line lysates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab123291 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab123291** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 86 kDa.

图标

<b>功能</b>	Calmodulin-binding protein which may function as scaffolding or signaling protein and may play a role in dendritic Ca(2+) signaling.
<b>组织特异性</b>	Preferentially expressed in brain.
<b>序列相似性</b>	Belongs to the WD repeat striatin family. Contains 6 WD repeats.
<b>细胞定位</b>	Cytoplasm. Membrane. Cell projection, dendritic spine. CTTNBP2-binding may regulate dendritic spine distribution.

## 图片



Anti-Striatin antibody (ab123291) at 1/100 dilution + NCI-H292 cell line lysates at 35 µg

Developed using the ECL technique.

**Predicted band size:** 86 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors