


Product datasheet

Anti-Syntaxin 5A antibody ab96185

1 图像

概述

产品名称	Anti-Syntaxin 5A抗体
描述	兔多克隆抗体to Syntaxin 5A
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Cow 
免疫原	Recombinant fragment containing a sequence corresponding to a region within amino acids 1-223 of Human Syntaxin 5A (AAH02645).
阳性对照	A431 whole cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab96185** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 30 kDa.

## 靶标

功能	Mediates endoplasmic reticulum to Golgi transport.
序列相似性	Belongs to the syntaxin family. Contains 1 t-SNARE coiled-coil homology domain.
细胞定位	Endoplasmic reticulum-Golgi intermediate compartment membrane. Golgi apparatus membrane.

## 图片



Anti-Syntaxin 5A antibody (ab96185) at  
1/5000 dilution + A431 whole cell lysate at 30  
 $\mu$ g

**Predicted band size:** 30 kDa

12% SDS-PAGE

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors