

## Product datasheet

# Anti-TEX11 antibody ab99461

### 1 References

#### 概述

**产品名称** Anti-TEX11抗体  
**描述** 山羊多克隆抗体to TEX11

该产品是 **快线抗体**，通过亲和纯化得到。对抗原肽有高灵敏度及特异性 (ELISA测试)。

[Read the terms of use »](#)

**宿主** Goat  
**种属反应性** 与反应: Human  
 预测可用于: Cynomolgus monkey

**免疫原** Synthetic peptide: HNFLKQTGTFSDSC, corresponding to internal sequence amino acids 722-736 of Human TEX11 (NP\_001003811.1; NP\_112566.2).  
[Run BLAST with ExPASy](#) [Run BLAST with NCBI](#)

#### 性能

**形式** Liquid  
**存放说明** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.  
**存储溶液** Preservative: 0.02% Sodium Azide  
 Constituents: 0.5% BSA, Tris buffered saline, pH 7.3  
**纯度** Immunogen affinity purified  
**纯化说明** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.  
**克隆** 多克隆  
**同种型** IgG

#### 应用

Fast track antibodies constitute a diverse group of products that have been released to accelerate your research, but are not yet fully characterized. They have all been affinity purified and show high titre values against the immunizing peptide (by ELISA).

[Fast track terms of use](#)

应用	Ab 评论	说明
----	-------	----

WB Use at an assay dependent dilution. Preliminary experiments in Human testis lysates gave no specific signal but low background (at antibody concentration up to 2 µg/ml). We would welcome your feedback if you have tested this antibody in more suitable positive controls.

## 靶标

**相关性** This gene is X-linked and is expressed in only male germ cells. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene

**细胞定位** Nuclear, Chromosome

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors