

Anti-TEX11 antibody ab99461

2 References

概述

产品名称

Anti-TEX11抗体

描述

山羊多克隆抗体to TEX11

① 该产品是 **快线抗体**，通过亲和纯化得到。对抗原肽有高灵敏度及特异性(ELISA测试)。

[Read the terms of use »](#)

宿主

Goat

种属反应性

与反应: Human

预测可用于: Cynomolgus monkey 

免疫原

Synthetic peptide: HNFLKQTGTFSDSC, corresponding to internal sequence amino acids 722-736 of Human TEX11 (NP_001003811.1; NP_112566.2).

 [Run BLAST with](#)  [Expasy](#)  [Run BLAST with](#)  [NCBI](#)

性能

形式

Liquid

存放说明

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

存储溶液

pH: 7.30

Preservative: 0.02% Sodium azide

Constituents: 0.5% BSA, Tris buffered saline

纯度

Immunogen affinity purified

纯化说明

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

克隆

多克隆

同种型

IgG

应用

Fast track antibodies constitute a diverse group of products that have been released to accelerate your research, but are not yet fully characterized. They have all been affinity purified and show high titre values against the immunizing peptide (by ELISA).

[Fast track terms of use](#)

应用	Ab评论	说明
WB		Use at an assay dependent concentration. A 1 hour primary incubation is recommended for this product. Preliminary experiments in Human testis lysates gave no specific signal but low background (at antibody concentration up to 2 µg/ml). We would welcome your feedback if you have tested this antibody in more suitable positive controls.

靶标

相关性 This gene is X-linked and is expressed in only male germ cells. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene

细胞定位 Nuclear, Chromosome

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors