

## Product datasheet

### Anti-TL1A antibody ab37169

#### 1 图像

#### 概述

产品名称	Anti-TL1A抗体
描述	兔多克隆抗体to TL1A
特异性	Does not crossreact with TL1.
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Human
免疫原	Recombinant full length TL1A (Human)
阳性对照	Mouse lung cell lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: PBS
纯度	IgG fraction
纯化说明	DEAE purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab37169** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 0.5 - 1 µg/ml. Detects a band of approximately 28 kDa (predicted molecular weight: 28 kDa).

#### 靶标

<b>功能</b>	Receptor for TNFRSF25 and TNFRSF6B. Mediates activation of NF-kappa-B. Inhibits vascular endothelial growth and angiogenesis (in vitro). Promotes activation of caspases and apoptosis.
<b>组织特异性</b>	Specifically expressed in endothelial cells. Detected in monocytes, placenta, lung, liver, kidney, skeletal muscle, pancreas, spleen, prostate, small intestine and colon.
<b>序列相似性</b>	Belongs to the tumor necrosis factor family.
<b>细胞定位</b>	Secreted and Membrane.

## 图片



**Predicted band size : 28 kDa**

Western Blot of human placenta tissue labeling TL1A with Anti-TL1A antibody (ab37169) at 2µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors