

Product datasheet

Anti-TMEM25 antibody ab121362

3 图像

概述

产品名称	Anti-TMEM25抗体
描述	兔多克隆抗体to TMEM25
宿主	Rabbit
经测试应用	适用于: ICC/IF, IHC-P, WB
种属反应性	与反应: Human
免疫原	antigen, corresponding to amino acids 262-364 of Human TMEM25 (Uniprot: Q86YD3)
阳性对照	Human Liver tissue and Liver lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab121362** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.
IHC-P		1/10 - 1/20. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

应用	Ab评论	说明
----	------	----


WB 1/250 - 1/500.

靶标

序列相似性 Contains 1 Ig-like (immunoglobulin-like) domain.

细胞定位 Secreted and Cell membrane.

图片



All lanes : Anti-TMEM25 antibody (ab121362) at 1/250 dilution

Lane 1 : RT-4 lysate


Lane 2 : U-251 MG lysate

Lane 3 : Human Plasma lysate

Lane 4 : Liver lysate


Lane 5 : Tonsil lysate

Western blot - Anti-TMEM25 antibody (ab121362)



ab121362, at 1/10 dilution, staining Paraffin-embedded Human Liver tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TMEM25 antibody (ab121362)



Immunofluorescent staining of human cell line A-431 shows positivity in vesicles.

Immunocytochemistry - Anti-TMEM25 antibody
(ab121362)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors