

Anti-TNFRSF14/HVEM antibody [MM0332-5S1] ab89479

2 References

概述

产品名称	Anti-TNFRSF14/HVEM抗体[MM0332-5S1]
描述	小鼠单克隆抗体[MM0332-5S1] to TNFRSF14/HVEM
宿主	Mouse
种属反应性	与反应: Human
免疫原	Recombinant fragment corresponding to Human TNFRSF14/HVEM (extracellular).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituent: PBS
纯度	Protein G purified
纯化说明	Purified from tissue culture supernatant and 0.2 µm filter sterilised.
克隆	单克隆
克隆编号	MM0332-5S1
同种型	IgG1

靶标

功能	Receptor for BTLA. Receptor for TNFSF14/LIGHT and homotrimeric TNFSF1/lymphotoxin-alpha. Involved in lymphocyte activation. Plays an important role in HSV pathogenesis because it enhanced the entry of several wild-type HSV strains of both serotypes into CHO cells, and mediated HSV entry into activated human T-cells.
----	---

组织特异性	Widely expressed, with the highest expression in lung, spleen and thymus.
序列相似性	Contains 3 TNFR-Cys repeats.
翻译后修饰	N-glycosylated.
细胞定位	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors