


Product datasheet

Anti-TNNI3K antibody ab96684

2 图像

概述

产品名称	Anti-TNNI3K抗体
描述	兔多克隆抗体to TNNI3K
经测试应用	适用于: WB, ICC/IF
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Cow 
免疫原	Recombinant fragment containing a sequence corresponding to a region within amino acids 187-413 of Human TNNI3K (NP_057062).
阳性对照	293T; A431; H1299; HeLaS3; HepG2 or Raji cell lysates

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab96684** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 93 kDa.
ICC/IF		1/100 - 1/200.

靶标


<b>功能</b>	May play a role in cardiac physiology.
<b>组织特异性</b>	Highly expressed in both adult and fetal heart.
<b>序列相似性</b>	Belongs to the protein kinase superfamily. TKL Ser/Thr protein kinase family. MAP kinase kinase kinase subfamily. Contains 10 ANK repeats. Contains 1 protein kinase domain.
<b>翻译后修饰</b>	Autophosphorylated.
<b>细胞定位</b>	Nucleus. Cytoplasm. Expressed at lower levels in the cytoplasm.

### Anti-TNNI3K antibody 图像



Anti-TNNI3K antibody (ab96684) at 1/5000 dilution + HeLa whole cell lysate at 30 µg

**Predicted band size : 93 kDa**



ab96684, at 1/200 dilution, staining TNNI3K in paraformaldehyde-fixed HeLa cells by Immunofluorescence. Lower image is merged with DNA probe.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors