

## Product datasheet

# TOM1L2 peptide (462-475) ab45647

### 概述

**产品名称** TOM1L2肽(462-475)

### 描述

**性质** Synthetic

### 氨基酸序列

**种属** Human

**序列** C-RKKPERSEDALFAL

**氨基酸** 462 to 475

### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab45647** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**应用** Blocking - Blocking peptide for Anti-TOM1L2 antibody ([ab24369](#))

**形式** Liquid

### 制备和贮存

**稳定性和存储** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### 常规信息

**功能** Probable role in protein transport. May regulate growth factor-induced mitogenic signaling.

**组织特异性** Ubiquitously expressed with higher expression in heart and skeletal muscle.

**序列相似性** Belongs to the TOM1 family.  
Contains 1 GAT domain.  
Contains 1 VHS domain.

**结构域** The GAT domain mediates interaction with TOLLIP.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors