

Product datasheet

Anti-TWEAKR antibody [ITEM-4] (Phycoerythrin) ab93555

2 图像

概述

产品名称	Anti-TWEAKR抗体[ITEM-4] (Phycoerythrin)
描述	小鼠单克隆抗体[ITEM-4] to TWEAKR (Phycoerythrin)
宿主	Mouse
偶联物	Phycoerythrin. Ex: 488nm, Em: 575nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse, Human
免疫原	Human TWEAKR transfected P815.
阳性对照	Human TWEAKR transfected L5178Y cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS, 150mM Sodium chloride, pH 7.2
纯度	Protein G purified
克隆	单克隆
克隆编号	ITEM-4
同种型	IgG2b
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab93555** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

Flow Cyt

应用说明

Flow Cyt: Use at a concentration of 0.5 - 1 µg/test.

A test is defined as the amount (μg) of antibody that will stain a cell sample in a final volume of $100\ \mu\text{L}$. Cell number should be determined empirically but can range from 10^5 to 10^8 cells/test.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能

Receptor for TNFSF12/TWEAK. Weak inducer of apoptosis in some cell types. Promotes angiogenesis and the proliferation of endothelial cells. May modulate cellular adhesion to matrix proteins.

组织特异性

Highly expressed in heart, placenta and kidney. Intermediate expression in lung, skeletal muscle and pancreas.


序列相似性

Contains 1 TNFR-Cys repeat.

细胞定位


Membrane.

图片



ab93555 at $0.5\ \mu\text{g}$ dilution staining TWEAKR in non-transfected L5178Y cells by Flow Cytometry (shown in blue).

Flow Cytometry - Anti-TWEAKR antibody [ITEM-4]
(Phycoerythrin) (ab93555)



ab93555 at $0.5\ \mu\text{g}$ dilution staining TWEAKR in Human TWEAKR transfected L5178Y cells by Flow Cytometry (shown in blue).

Flow Cytometry - Anti-TWEAKR antibody [ITEM-4]
(Phycoerythrin) (ab93555)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors