


Product datasheet

Anti-Tetraspanin 9 antibody ab113775

2 图像

概述

产品名称	Anti-Tetraspanin 9抗体
描述	兔多克隆抗体to Tetraspanin 9
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Mouse, Human 预测可用于: Rat 
免疫原	KLH conjugated synthetic peptide of an internal region of Human Tetraspanin 9.
阳性对照	EL4 cell lysate, Human Thyroid tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab113775** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 27 kDa.
IHC-P		Use a concentration of 5 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

靶标

组织特异性	Expressed in megakaryocytes and platelets (at protein level).
序列相似性	Belongs to the tetraspanin (TM4SF) family.
翻译后修饰	Glycosylated.
细胞定位	Membrane. Colocalizes with GP6 in tetraspanin microdomains on the platelet surface.

图片




All lanes : Anti-Tetraspanin 9 antibody (ab113775) at 1 μ g/ml

Lane 1 : EL4 cell lysate

Lane 2 : EL4 cell lysate with blocking peptide

Predicted band size : 27 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Tetraspanin 9 antibody (ab113775)

Staining of Human Tetraspanin 9 in a Formalin/PFA-fixed paraffin-embedded section of Human Thyroid using ab113775 at a concentration of 5 μ g/ml.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors