


Product datasheet

Anti-Tetraspanin 9 antibody ab113775

2 图像

概述

产品名称	Anti-Tetraspanin 9抗体
描述	兔多克隆抗体to Tetraspanin 9
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Mouse, Human 预测可用于: Rat 
免疫原	KLH conjugated synthetic peptide of an internal region of Human Tetraspanin 9.
阳性对照	EL4 cell lysate, Human Thyroid tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab113775** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 27 kDa.

应用	Ab评论	说明
IHC-P		Use a concentration of 5 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## 靶标

组织特异性	Expressed in megakaryocytes and platelets (at protein level).
序列相似性	Belongs to the tetraspanin (TM4SF) family.
翻译后修饰	Glycosylated.
细胞定位	Membrane. Colocalizes with GP6 in tetraspanin microdomains on the platelet surface.

## 图片




Western blot - Anti-Tetraspanin 9 antibody (ab113775)

**All lanes :** Anti-Tetraspanin 9 antibody (ab113775) at 1 µg/ml

**Lane 1 :** EL4 cell lysate

**Lane 2 :** EL4 cell lysate with blocking peptide

**Predicted band size:** 27 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Tetraspanin 9 antibody (ab113775)

Staining of Human Tetraspanin 9 in a Formalin/PFA-fixed paraffin-embedded section of Human Thyroid using ab113775 at a concentration of 5 µg/ml.

---

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors