

Product datasheet

Anti-UGT2B11 antibody ab116019

1 图像

概述

|       |  |
|-------|--|
| 产品名称  | Anti-UGT2B11抗体   |
| 描述    | 兔多克隆抗体to UGT2B11   |
| 经测试应用 | 适用于: WB  |
| 种属反应性 | 与反应: Human   |
| 免疫原   | Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 428-458 of Human UGT2B11 (NP_001064.1). |
| 阳性对照  | NCI-H460 cell line lysates.  |

性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.                    |
| 存储溶液 | Preservative: 0.09% Sodium azide<br>Constituent: 99% PBS                                    |
| 纯度   | Immunogen affinity purified   |
| 纯化说明 | ab116019 is purified through a protein A column, followed by peptide affinity purification. |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

应用

Our [Abpromise guarantee](#) covers the use of **ab116019** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明   |
|----|------|--|
| WB |      | 1/100 - 1/500. Predicted molecular weight: 61 kDa. |

靶标

|              |   |
|--------------|---|
| <b>功能</b>    | UDPGT is of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds. |
| <b>组织特异性</b> | Widely expressed.   |
| <b>序列相似性</b> | Belongs to the UDP-glycosyltransferase family.  |
| <b>细胞定位</b>  | Microsome membrane. Endoplasmic reticulum membrane.   |

## 图片



Anti-UGT2B11 antibody (ab116019) at 1/100 dilution + NCI-H460 cell line lysates at 35 µg

**Predicted band size : 61 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors