

Product datasheet

Anti-Urocortin antibody ab80358

1 图像

概述

产品名称	Anti-Urocortin抗体
描述	鸡多克隆抗体to Urocortin
宿主	Chicken
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant fragment, corresponding to amino acids 68-124 of human Urocortin derived from <i>E. coli</i>
阳性对照	<i>E. coli</i> derived Urocortin fusion protein.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgY

应用

Our [Abpromise guarantee](#) covers the use of **ab80358** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/3000. Predicted molecular weight: 13 kDa.

靶标

功能	Acts in vitro to stimulate the secretion of adrenocorticotrophic hormone (ACTH). Binds with high affinity to CRF receptor types 1, 2-alpha, and 2-beta.
序列相似性	Belongs to the sauvagine/corticotropin-releasing factor/urotensin I family.
细胞定位	Secreted.

图片



Anti-Urocortin antibody (ab80358) at 1/5000 dilution + *E. coli* derived Urocortin fusion protein at 0.4 µg

Secondary

goat-anti-IgY/Fc-HRP at 1/10000 dilution

Developed using the ECL technique.

Predicted band size: 13 kDa

Observed band size: 13 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors