

Product datasheet

Anti-Urokinase antibody (FITC) ab28223

概述

产品名称	Anti-Urokinase抗体(FITC)
描述	兔多克隆抗体to Urokinase (FITC)
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Full length protein: uPA purified from human urine.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: None Constituents: 0.1M Sodium chloride, 0.05M Sodium phosphate, 1mM EDTA. pH 6.6.
纯化说明	< 98% pure by SDS PAGE.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab28223** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/1000. Predicted molecular weight: 48 kDa.
ELISA		Use at an assay dependent dilution.

靶标

功能	Specifically cleave the zymogen plasminogen to form the active enzyme plasmin.
----	--

组织特异性	Expressed in the prostate gland and prostate cancers.
序列相似性	Belongs to the peptidase S1 family. Contains 1 EGF-like domain. Contains 1 kringle domain. Contains 1 peptidase S1 domain.
翻译后修饰	Phosphorylation of Ser-158 and Ser-323 abolishes proadhesive ability but does not interfere with receptor binding.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors