

Product datasheet

Anti-Urokinase antibody (FITC) ab28223

概述

产品名称	Anti-Urokinase抗体(FITC)
描述	兔多克隆抗体to Urokinase (FITC)
宿主	Rabbit
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Full length protein: uPA purified from human urine.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: None Constituents: 0.1M Sodium chloride, 0.05M Sodium phosphate, 1mM EDTA. pH 6.6.
纯化说明	< 98% pure by SDS PAGE.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab28223** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/1000. Predicted molecular weight: 48 kDa.
ELISA		Use at an assay dependent dilution.

靶标

功能	Specifically cleave the zymogen plasminogen to form the active enzyme plasmin.
组织特异性	Expressed in the prostate gland and prostate cancers.
序列相似性	Belongs to the peptidase S1 family. Contains 1 EGF-like domain. Contains 1 kringle domain. Contains 1 peptidase S1 domain.
翻译后修饰	Phosphorylation of Ser-158 and Ser-323 abolishes proadhesive ability but does not interfere with receptor binding.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors