

## Product datasheet

# Anti-Vibrio parahaemolyticus K24 antibody ab79501

### 概述

产品名称	Anti-Vibrio parahaemolyticus K24抗体
描述	兔多克隆抗体to Vibrio parahaemolyticus K24
特异性	Specific to <i>Vibrio parahaemolyticus</i> K24 antigen.
经测试应用	适用于: Agglutination
种属反应性	
免疫原	Tissue/ cell preparation of <i>Vibrio parahaemolyticus</i> K24: (heat inactivated)
阳性对照	Antigenic suspension
常规说明	Reagents should be allowed to stand at 15°C – 25°C for at least 30 minutes prior to use.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab79501** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

Agglutination

应用说明	Agg: See attached protocol.  Not yet tested in other applications.
------	--

Optimal dilutions/concentrations should be determined by the end user.

## 靶标

---

### 相关性

*Vibrio parahaemolyticus* is an oxidase positive, facultatively aerobic, non spore forming, Gram negative bacillus with a single polar flagellum. When ingested it causes gastrointestinal illness in humans, with ingestion via raw or undercooked seafood (typically oysters) being the main route of infection.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors