

Product datasheet

Anti-Vibrio parahaemolyticus K24 antibody ab79501

概述

产品名称	Anti-Vibrio parahaemolyticus K24抗体
描述	兔多克隆抗体to Vibrio parahaemolyticus K24
宿主	Rabbit
特异性	Specific to <i>Vibrio parahaemolyticus</i> K24 antigen.
经测试应用	适用于: Agglutination
种属反应性	
免疫原	Tissue/ cell preparation of <i>Vibrio parahaemolyticus</i> K24: (heat inactivated)
阳性对照	Antigenic suspension
常规说明	Reagents should be allowed to stand at 15°C – 25°C for at least 30 minutes prior to use.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab79501** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Agglutination		Use at an assay dependent concentration. See attached protocol.

靶标

相关性

Vibrio parahaemolyticus is an oxidase positive, facultatively aerobic, non spore forming, Gram negative bacillus with a single polar flagellum. When ingested it causes gastrointestinal illness in humans, with ingestion via raw or undercooked seafood (typically oysters) being the main route of infection.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors