

Product datasheet

Anti-Vibrio parahaemolyticus K55 antibody ab78959

概述

产品名称	Anti-Vibrio parahaemolyticus K55抗体
描述	兔多克隆抗体to Vibrio parahaemolyticus K55
特异性	ab78959 is specific to Vibrio parahaemolyticus K55 antigen
经测试应用	适用于: Agglutination
种属反应性	Reacts with Vibrio parahaemolyticus K55
免疫原	Tissue/ cell preparation of Vibrio parahaemolyticus K55 (heat inactivated)
阳性对照	Antigenic suspension.
常规说明	Reagents should be allowed to stand at 15°C - 25°C for at least 30 minutes prior to use.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.08% Sodium Azide Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab78959** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Agglutination		Use at an assay dependent concentration. See attached protocol.

相关性

Vibrio parahaemolyticus, a member of the *Vibrio* genus, is a curved, rod shaped, Gram negative bacterium, found in brackish saltwater, which, when ingested, causes gastrointestinal illness in humans. *Vibrio parahaemolyticus* is oxidase positive, facultatively aerobic, non spore forming and is motile, with a single, polar flagellum. Most people become infected with *Vibrio parahaemolyticus* by eating raw or undercooked shellfish, particularly oysters. Less commonly, this organism can cause an infection in the skin when an open wound is exposed to warm seawater.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors