

Product datasheet

Anti-WIF1 antibody [RM0144-3M51] ab86496

概述

| | |
|-------|------------------------------|
| 产品名称 | Anti-WIF1抗体[RM0144-3M51] |
| 描述 | 大鼠单克隆抗体[RM0144-3M51] to WIF1 |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Mouse |
| 免疫原 | Recombinant mouse WIF1 |

性能

| | |
|------|--|
| 形式 | Lyophilised:Reconstitute the antibody with 200 µl sterile PBS for a final antibody concentration of 0.5 mg/ml. |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None Constituents: PBS |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | IgG fraction of culture supernatant purified by Protein A/G affinity chromatography and 0.2 µm filtered. |
| 克隆 | 单克隆 |
| 克隆编号 | RM0144-3M51 |
| 同种型 | IgG2 |

应用

Our [Abpromise guarantee](#) covers the use of **ab86496** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|---|
| WB | | 1/500 - 1/1000. Predicted molecular weight: 42 kDa. |

靶标

| | |
|-------|--|
| 功能 | Binds to WNT proteins and inhibits their activities. May be involved in mesoderm segmentation. |
| 序列相似性 | Contains 5 EGF-like domains. Contains 1 WIF domain. |
| 细胞定位 | Secreted. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors