

## Product datasheet

# Anti-WIF1 antibody [RM0144-3M51] ab86496

### 概述

产品名称	Anti-WIF1抗体[RM0144-3M51]
描述	大鼠单克隆抗体[RM0144-3M51] to WIF1
经测试应用	适用于: WB
种属反应性	与反应: Mouse
免疫原	Recombinant mouse WIF1

### 性能

形式	Lyophilised:Reconstitute the antibody with 200 µl sterile PBS for a final antibody concentration of 0.5 mg/ml.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	IgG fraction of culture supernatant purified by Protein A/G affinity chromatography and 0.2 µm filtered.
克隆	单克隆
克隆编号	RM0144-3M51
同种型	IgG2

### 应用

Our [Abpromise guarantee](#) covers the use of **ab86496** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 42 kDa.

### 靶标

<b>功能</b>	Binds to WNT proteins and inhibits their activities. May be involved in mesoderm segmentation.
<b>序列相似性</b>	Contains 5 EGF-like domains. Contains 1 WIF domain.
<b>细胞定位</b>	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors