

## Product datasheet

# Anti-XIRP1 antibody ab116138

★☆☆☆☆ 1 Abreviews 1 图像

### 概述

产品名称	Anti-XIRP1抗体
描述	兔多克隆抗体to XIRP1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 1337-1367 of Human XIRP1 (NP_919269.2).
阳性对照	NCI-H460 cell line lysates.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab116138 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab116138** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB	★☆☆☆☆	1/100 - 1/500. Predicted molecular weight: 199 kDa.

### 标签

<b>功能</b>	Protects actin filaments from depolymerization.
<b>组织特异性</b>	Isoform A, isoform B and isoform C are expressed in heart.
<b>序列相似性</b>	Belongs to the Xin family. Contains 17 Xin repeats.
<b>结构域</b>	Xin repeats bind F-actin.
<b>细胞定位</b>	Cell junction. Colocalizes with actin stress fibers.

## 图片



Anti-XIRP1 antibody (ab116138) at 1/100 dilution + NCI-H460 cell line lysate at 35 µg

**Predicted band size:** 199 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors