

Product datasheet

Anti-YIPF3 antibody ab121456

3 图像

概述

产品名称	Anti-YIPF3抗体
描述	兔多克隆抗体to YIPF3
宿主	Rabbit
经测试应用	适用于: WB, IHC-P, ICC/IF
种属反应性	与反应: Human
免疫原	antigen sequence, corresponding to amino acids 302-350 (GILDTLEGPN IPPIQRVPRD IPAMLPAARL PTTVLNATAK AVAVTLQSH) of Human YIPF3. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	Transfected HEK293T Over-expression Lysate; Human colon tissue; A431 cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab121456** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/250 - 1/500.

应用	Ab评论	说明
IHC-P		1/10 - 1/20. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		Use a concentration of 1 - 4 µg/ml. Fixation/Permeabilization: PFA/Triton X-100

## 靶标

功能	May play a role in hematopoiesis.
组织特异性	Expressed by nucleated hematopoietic cells (at protein level).
序列相似性	Belongs to the YIP1 family.
发展阶段	Expressed in fetal liver (at protein level). Expressed in embryonic hematopoiesis sites.
翻译后修饰	Not glycosylated.
细胞定位	Cell membrane. Cytoplasm. Localization to the cytoplasm or to the cell membrane is developmentally and ontogenetically regulated.

## 图片




**All lanes** : Anti-YIPF3 antibody (ab121456) at 1/250 dilution

**Lane 1** : Negative control (vector only transfected HEK293T lysate)


**Lane 2** : YIPF3-transfected HEK293T overexpression Lysate (co-expressed with a C terminal myc-DDK tag (~3.1 kDa))

Developed using the ECL technique.



ab121456, at a 1/15 dilution, staining YIPF3 in paraffin embedded Human colon tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-YIPF3 antibody (ab121456)



ab121456, at 4 µg/ml, staining YIPF3 (green) in Human A431 cells by Immunofluorescence.

Immunocytochemistry/ Immunofluorescence - Anti-YIPF3 antibody (ab121456)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors