

Product datasheet

Anti-YTHDF2 antibody ab88809

2 图像

概述

产品名称	Anti-YTHDF2抗体
描述	小鼠多克隆抗体to YTHDF2
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length Human YTHDF2 (NP_057342.2)
阳性对照	Human liver lysate. YTHDF2 transfected 293T cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Constituent: PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab88809** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

应用说明 WB: 1/500 - 1/1000. Predicted molecular weight: Isoform 1, 62 kDa; Isoform 2, 57 kDa.

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

靶标

序列相似性

Contains 1 YTH domain.

图片




Western blot - YTHDF2 antibody (ab88809)

Anti-YTHDF2 antibody (ab88809) at 1/500 dilution + Human liver lysate at 50 μ g

Predicted band size : 62 kDa

Observed band size : 49 kDa



Western blot - YTHDF2 antibody (ab88809)

All lanes : Anti-YTHDF2 antibody (ab88809) at 1/500 dilution

Lane 1 : YTHDF2 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 μ g per lane.

Predicted band size : 62 kDa

Observed band size : 70 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors