

Product datasheet

Anti-Z DNA binding protein antibody ab32927

1 图像

概述

产品名称	Anti-Z DNA binding蛋白抗体
描述	兔多克隆抗体to Z DNA binding蛋白
宿主	Rabbit
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Human
免疫原	Synthetic peptide within Human Z DNA binding protein aa 144-193 (internal sequence). The exact sequence is proprietary. Sequence: LYRMKSRHLLDMDEQSKAWTIYRPEDSGRRAKSASIIYQHNPINMICQNG Database link: Q9H171
阳性对照	Jurkat cell lysate.

 [Run BLAST with](#)

 [Run BLAST with](#)

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab32927** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
ELISA		1/312500.
WB		Use a concentration of 5 µg/ml. Detects a band of approximately 46 kDa (predicted molecular weight: 46 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能	May play a role in host defense against tumors and pathogens. Binds Z-DNA.
组织特异性	Highly expressed in lymphatic tissues including lymph node, leukocytes, tonsil, bone marrow and spleen. Expressed to a lesser extent in thymus, lung and liver.
序列相似性	Contains 2 DRADA repeats.

图片



Anti-Z DNA binding protein antibody (ab32927) at 5 µg/ml + Jurkat cell lysate

Secondary

HRP conjugated anti-Rabbit IgG
1: 50,000 - 100,000

Predicted band size: 46 kDa

Observed band size: 46 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors