

Product datasheet

Anti-Z DNA binding protein antibody ab32927

1 图像

概述

<b>产品名称</b>	Anti-Z DNA binding protein抗体
<b>描述</b>	兔多克隆抗体to Z DNA binding protein
<b>经测试应用</b>	适用于: ELISA, WB
<b>种属反应性</b>	与反应: Human
<b>免疫原</b>	Synthetic peptide corresponding to Human Z DNA binding protein aa 144-193 (internal sequence). Sequence: LYRMKSRHLLDMDEQSKAWTIYRPEDSGRRAKSASIIYQHNPINMIC  Database link: <a href="#">Q9H171</a>
	<a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
<b>阳性对照</b>	Jurkat cell lysate.

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>存储溶液</b>	Preservative: None Constituents: 2% Sucrose, PBS.
<b>纯度</b>	Protein A purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab32927** in the following tested applications.  
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
ELISA		1/312500.
WB		Use a concentration of 5 µg/ml. Detects a band of approximately 46 kDa (predicted molecular weight: 46 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

<b>功能</b>	May play a role in host defense against tumors and pathogens. Binds Z-DNA.
<b>组织特异性</b>	Highly expressed in lymphatic tissues including lymph node, leukocytes, tonsil, bone marrow and spleen. Expressed to a lesser extent in thymus, lung and liver.
<b>序列相似性</b>	Contains 2 DRADA repeats.

## 图片



Anti-Z DNA binding protein antibody (ab32927) at 5 µg/ml + Jurkat cell lysate

### Secondary

HRP conjugated anti-Rabbit IgG 1: 50,000 - 100,000

**Predicted band size** : 46 kDa

**Observed band size** : 46 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors