



Product datasheet

Anti-Z DNA binding protein antibody ab32927

1 图像

概述

| | |
|--------------|---|
| 产品名称 | Anti-Z DNA binding蛋白抗体 |
| 描述 | 兔多克隆抗体to Z DNA binding蛋白 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: ELISA, WB |
| 种属反应性 | 与反应: Human |
| 免疫原 | Synthetic peptide within Human Z DNA binding protein aa 144-193 (internal sequence). The exact sequence is proprietary. Sequence: LYRMKSRHLLDMDEQSKAWTIYRPEDSGRRAKSASIIYQHNPINMICQNG Database link: Q9H171 |
| |  Run BLAST with  Run BLAST with |
| 阳性对照 | Jurkat cell lysate. |

性能

| | |
|-------------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Preservative: None Constituents: 2% Sucrose, PBS. |
| 纯度 | Protein A purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab32927** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab 评论 | 说明 |
|-------|-------|---|
| ELISA | | 1/312500. |
| WB | | Use a concentration of 5 µg/ml. Detects a band of approximately 46 kDa (predicted molecular weight: 46 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

靶标

| | |
|--------------|---|
| 功能 | May play a role in host defense against tumors and pathogens. Binds Z-DNA. |
| 组织特异性 | Highly expressed in lymphatic tissues including lymph node, leukocytes, tonsil, bone marrow and spleen. Expressed to a lesser extent in thymus, lung and liver. |
| 序列相似性 | Contains 2 DRADA repeats. |

图片



Anti-Z DNA binding protein antibody (ab32927) at 5 µg/ml + Jurkat cell lysate

Secondary

HRP conjugated anti-Rabbit IgG
1: 50,000 - 100,000

Predicted band size: 46 kDa

Observed band size: 46 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors