






Anti-ZHX2 antibody ab66595

1 图像


概述		
产品名称	Anti-ZHX2抗体	
描述	兔多克隆抗体to ZHX2	
宿主	Rabbit	
经测试应用	适用于: WB	
种属反应性	与反应: Mouse	
	预测可用于: Rat, Rabbit, Cow, Cat, Dog, Pig 	
免疫原	A region within synthetic peptide: GMDFVEVTV GEEDAISEKW GSWSRRVAEG TVERADSDSD STPAEAGQA, corresponding to C terminal amino acids 788-836 of Mouse ZHX2	
	 Run BLAST with   Run BLAST with 	
阳性对照	Sp2/0 lysate	
性能		
形式	Liquid	
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.	
存储溶液	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS	
纯度	Protein A purified	
克隆	多克隆	
同种型	IgG	
应用		
The Abpromise guarantee		Abpromise™ 承诺保证使用ab66595于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。		
应用	Ab评论	说明

应用	Ab评论	说明
WB		Use a concentration of 2.5 µg/ml. Detects a band of approximately 90 kDa (predicted molecular weight: 92 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能	Acts as a transcriptional repressor.
组织特异性	Ubiquitously expressed.
序列相似性	Belongs to the ZHX family. Contains 2 C2H2-type zinc fingers. Contains 4 homeobox DNA-binding domains.
细胞定位	Nucleus.

图片



Anti-ZHX2 antibody (ab66595) at 2.5 µg/ml + SP 2/0 lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 92 kDa

Observed band size: 90 kDa

Gel concentration: 12%.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors