

Product datasheet

Anti-ZNF280D antibody ab114972

2 图像

概述

产品名称	Anti-ZNF280D抗体
描述	兔多克隆抗体to ZNF280D
宿主	Rabbit
经测试应用	适用于: WB, IP
种属反应性	与反应: Human 预测可用于: Horse, Dog, Turkey, Chimpanzee, Rhesus monkey, Gorilla, Opossum, Orangutan, Bat ▲
免疫原	Synthetic peptide, corresponding to a region between residue 879 and 929 of Human ZNF280D (NP_060131.2).
阳性对照	HeLa, 293T and Jurkat whole cell lysates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% Tris citrate/phosphate Note: pH 7 to 8
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab114972** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/1000 - 1/5000. Predicted molecular weight: 109 kDa.

应用	Ab评论	说明
IP		Use at 2-10 µg/mg of lysate.

## 靶标

功能	May function as a transcription factor.
序列相似性	Contains 5 C2H2-type zinc fingers.
细胞定位	Nucleus.

## 图片



**All lanes :** Anti-ZNF280D antibody (ab114972) at 0.4 µg/ml

**Lane 1 :** HeLa whole cell lysate at 50 µg

**Lane 2 :** HeLa whole cell lysate at 15 µg


**Lane 3 :** 293T whole cell lysate at 50 µg

**Lane 4 :** Jurkat whole cell lysate at 50 µg

Developed using the ECL technique.

**Predicted band size:** 109 kDa

**Exposure time:** 30 seconds



Detection of ZNF280D in Immunoprecipitates of HeLa whole cell lysate (1 mg for IP, 20% of IP loaded) using ab114972 at 6 µg/mg lysate for IP and at 1 µg/ml for subsequent Western blot detection. Detection: Chemiluminescence with an exposure time of 10 seconds.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors