

Product datasheet

Anti-ZNF398 antibody ab96651

2 图像

概述

产品名称	Anti-ZNF398抗体
描述	兔多克隆抗体to ZNF398
宿主	Rabbit
经测试应用	适用于: WB, ICC/IF
种属反应性	与反应: Human 预测可用于: Mouse, Rat
免疫原	Recombinant protein fragment corresponding to a region within amino acids 358 and 630 of ZNF398 (NP_733787).
阳性对照	293T, A431, H1299, HeLa, MOLT4 and Raji cell lines

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab96651** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 71 kDa.
ICC/IF		1/100 - 1/200.

靶标

功能	Function as a transcriptional activator.
序列相似性	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 9 C2H2-type zinc fingers. Contains 1 KRAB domain.
细胞定位	Nucleus.

图片




All lanes : Anti-ZNF398 antibody (ab96651)
at 1/1000 dilution

Lane 1 : 293T whole cell lysate
Lane 2 : H1299 whole cell lysate

Lysates/proteins at 30 μ g per lane.

Predicted band size: 71 kDa

7.5% SDS PAGE



ab96651, at a 1/200 dilution, staining ZNF398
in paraformaldehyde fixed HeLa by
Immunofluorescence analysis. Lower image:
merged with DNA probe

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors