

Product datasheet

Anti-ZNF563 antibody ab70905

2 图像

概述

产品名称	Anti-ZNF563抗体
描述	小鼠多克隆抗体to ZNF563
经测试应用	适用于: ICC/IF, WB
种属反应性	与反应: Human
免疫原	Full length human ZNF563 protein (NP_660319.1).
阳性对照	ZNF563 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab70905** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		Use a concentration of 10 µg/ml.
WB		1/500 - 1/1000. Detects a band of approximately 48 kDa (predicted molecular weight: 55 kDa).

靶标

<b>功能</b>	May be involved in transcriptional regulation.
<b>序列相似性</b>	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 12 C2H2-type zinc fingers. Contains 1 KRAB domain.
<b>细胞定位</b>	Nucleus.

## 图片



**All lanes** : Anti-ZNF563 antibody (ab70905)  
at 1/500 dilution

**Lane 1** : ZNF563 transfected 293T cell lysate

**Lane 2** : Non transfected 293T cell lysate


Lysates/proteins at 25 µg per lane.

### Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at  
1/2500 dilution

**Predicted band size** : 55 kDa

**Observed band size** : 48 kDa



ab70905 at 10µg/ml staining ZNF563 in  
human HeLa cells by ICC/IF.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors