


Product datasheet

Anti-ZNF690 antibody ab117831

1 图像

概述

产品名称	Anti-ZNF690抗体
描述	兔多克隆抗体to ZNF690
经测试应用	适用于: IP 不适用于: WB
种属反应性	与反应: Human 预测可用于: Chimpanzee, Baboon, Rhesus monkey, Gorilla, Orangutan 
免疫原	Synthetic peptide, corresponding to a region within amino acids 600-650 of Human ZNF690 (NP_689668.3).
阳性对照	HeLa whole cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituent: 99% Tris citrate/phosphate Note: pH: 7-8
纯度	Immunogen affinity purified
纯化说明	ab117831 was affinity purified using an epitope specific to ZNF690 immobilized on solid support.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab117831** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IP		Use at 2-10 µg/mg of lysate.

应用说明

Is unsuitable for WB.

靶标

功能

May be involved in transcriptional regulation.

序列相似性

Belongs to the krueppel C2H2-type zinc-finger protein family.


Contains 6 C2H2-type zinc fingers.

Contains 1 SCAN box domain.

细胞定位

Nucleus.

图片



Immunoprecipitation - Anti-ZNF690 antibody
(ab117831)

Detection of Human ZNF690 in HeLa whole cell lysate (1 mg for IP, 20% of IP loaded) using ab117831 at 6 μ g/mg of lysate.

Subsequent WB detection was done using an other anti Human ZNF690 antibody.

Detection: Chemiluminescence with an exposure time of 3 minutes.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors