

Product datasheet

Anti-Rotavirus group specific antigen antibody [3C10] ab10012

概述

产品名称	Anti-Rotavirus group specific antigen抗体[3C10]
描述	小鼠单克隆抗体[3C10] to Rotavirus group specific antigen
宿主	Mouse
特异性	Rotavirus P42 inner-capsid antigen specific. Cross-reacts with monkey rotavirus (SA-11), porcine rotavirus (PP) and human rotaviruses (field strains).
经测试应用	适用于: WB
免疫原	Full length native protein (purified bovine rotavirus, cell culture adapted strain MR).
常规说明	Concentration varies from lot to lot and can be provided on request.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.4
纯度	Protein A purified
纯化说明	Purity tested by elcetrophoresis.
克隆	单克隆
克隆编号	3C10
骨髓瘤	Sp2/0
同种型	IgG2a

应用

Our [Abpromise guarantee](#) covers the use of **ab10012** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors