

Product datasheet

Anti-ZNF514 antibody ab101511

1 图像

概述

产品名称	Anti-ZNF514抗体
描述	兔多克隆抗体to ZNF514
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant fragment, corresponding to a region within amino acids 71-215 of Human ZNF514 (AAH74900).
阳性对照	HepG2 whole cell lysate; A549, HeLa and HCT116 cell lysates

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab101511** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 46 kDa.

靶标

功能	May be involved in transcriptional regulation.
----	--

序列相似性

Belongs to the krueppel C2H2-type zinc-finger protein family.


Contains 7 C2H2-type zinc fingers.

Contains 1 KRAB domain.

细胞定位

Nucleus.

图片



Anti-ZNF514 antibody (ab101511) at 1/2000
dilution + HepG2 whole cell lysate at 30 µg

Predicted band size : 46 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors