

## Product datasheet

# Anti-PR3 antibody [1F10] ab10198

### 概述

产品名称	Anti-PR3抗体[1F10]
描述	小鼠单克隆抗体[1F10] to PR3
宿主	Mouse
特异性	<a href="#">ab10197</a> , 10198, 10199, 10200 and 10201 detect different determinants of proteinase PR3 which have not been fully determined.
经测试应用	适用于: ELISA
种属反应性	与反应: Human
免疫原	Human proteinase PR3 purified leukocyte mass
常规说明	Concentration varies from lot to lot and can be provided on request.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.4
纯度	Protein A purified
纯化说明	Purified by chromatography on protein A sepharose. Purity tested by electrophoresis.
克隆	单克隆
克隆编号	1F10
骨髓瘤	Sp2/0
同种型	IgG2b

### 应用

Our [Abpromise guarantee](#) covers the use of **ab10198** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.

## 靶标

---

功能	Polymorphonuclear leukocyte serine protease that degrades elastin, fibronectin, laminin, vitronectin, and collagen types I, III, and IV (in vitro) and causes emphysema when administered by tracheal insufflation to hamsters.
序列相似性	Belongs to the peptidase S1 family. Elastase subfamily. Contains 1 peptidase S1 domain.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors