

## Product datasheet

# Anti-FUBP3 antibody ab102576

### 2 图像

#### 概述

产品名称	Anti-FUBP3抗体
描述	兔多克隆抗体to FUBP3
宿主	Rabbit
经测试应用	适用于: WB, IHC-P, ELISA
种属反应性	与反应: Human 预测可用于: Mouse, Rat
免疫原	Recombinant protein fragment containing a sequence corresponding to a region within amino acids 14 and 153 of Human FUBP3 (AAH07874; UniProt ID: Q96I24 isoform 2).
阳性对照	WB: NT2D1, IMR32, U87-MG and MCF7 cell lysates IHC-P: Cal27 xenograft

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab102576** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/5000 - 1/20000. Predicted molecular weight: for isoform 2, 29 kDa.
IHC-P		1/100 - 1/500. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.

## 靶标

功能	May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.
组织特异性	Detected in a number of cell lines.
序列相似性	Contains 4 KH domains.
翻译后修饰	Phosphorylated upon DNA damage, probably by ATM or ATR.
细胞定位	Nucleus.


## 图片



Anti-FUBP3 antibody (ab102576) at 1/10000 dilution + MCF7 whole cell lysate at 30  $\mu$ g

**Predicted band size:** for isoform 2, 29 kDa

7.5% SDS PAGE



Immunohistochemical analysis of FUBP3 in paraffin embedded Cal27 xenograft, using ab102576 at a 1/100 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - FUBP3 antibody (ab102576)

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors