

Product datasheet

Anti-TC2N antibody ab102981

2 图像

概述

产品名称	Anti-TC2N抗体
描述	兔多克隆抗体to TC2N
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Mouse, Human
免疫原	KLH conjugated synthetic peptide corresponding to a region between N terminal amino acids 9 - 39 of Human TC2N.
阳性对照	WB: Mouse NIH-3T3 cell line lysates IHC-P: Human testis tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	ab102981 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab102981** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 55 kDa.

应用	Ab评论	说明
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

靶标


相关性	The specific function of TC2N is not yet known.
细胞定位	Nuclear

图片



Anti-TC2N antibody (ab102981) at 1/100 dilution + mouse NIH-3T3 cell line lysates at 35 μ g

Predicted band size: 55 kDa



ab102981, at a 1/50 dilution, staining TC2N in formalin-fixed and paraffin-embedded Human testis tissue, which was peroxidase-conjugated to the secondary antibody, followed by DAB staining.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors