

Product datasheet

Anti-Retinoid X Receptor beta antibody ab103027

1 图像

概述

产品名称	Anti-Retinoid X Receptor beta 抗体
描述	兔多克隆抗体 to Retinoid X Receptor beta
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 262-292 of Human Retinoid X Receptor beta.
阳性对照	MDA-MB231 cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	ab103027 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab103027** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 57 kDa.

靶标


## 相关性

Retinoid X receptors (RXRs) are members of the steroid/thyroid hormone receptor superfamily of nuclear receptor proteins which exert their effects by binding to specific DNA response elements, thus regulating gene expression in target cells. The RXR subfamily consists of at least three similar genes, RXR alpha, RXR beta and RXR gamma, all of which control transcription of target genes mediated by retinoids. RXR beta controls expression of many genes that respond to hormones and vitamins, including thyroid hormone, estrogen, retinoids and vitamin D. RXR beta controls a wide array of genes because of its ability to heterodimerize with other hormone receptors including the thyroid hormone receptor (THR), vitamin D receptor (VDR) and retinoic acid receptor (RAR).

## 细胞定位

Nuclear

## 图片



Anti-Retinoid X Receptor beta antibody (ab103027) at 1/100 dilution + MDA-MB231 cell lysate at 35 µg

**Predicted band size: 57 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors