

Product datasheet

Anti-HUS1B antibody ab104713

1 图像

概述

产品名称	Anti-HUS1B抗体
描述	兔多克隆抗体to HUS1B
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 14-63 (ELFIHVS ^G TVARLAKVCVLRVRPDSL ^C FGPAGSGGLHEARLWCEVR ^Q QAF) of Human HUS1B (NP_683762). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Human liver cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab104713** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 0.5 µg/ml. Predicted molecular weight: 31 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

相关性

This protein is most closely related to HUS1, a component of a cell cycle checkpoint protein complex involved in cell cycle arrest in response to DNA damage. This protein can interact with the checkpoint protein RAD1 but not with RAD9. Over-expression of this protein has been shown to induce cell death, which suggests a related but distinct role of this protein, as compared to the HUS1.

图片



Anti-HUS1B antibody (ab104713) at 0.5 µg/ml
+ Human liver cell lysate at 10 µg

Predicted band size: 31 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors