

## Product datasheet

# Anti-LTBR antibody ab10493

### 概述

产品名称	Anti-LTBR抗体
描述	山羊多克隆抗体to LTBR
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Human
免疫原	Recombinant human lymphotoxin betaR extracellular domain, expressed in mouse NSO cells.
常规说明	Endotoxin level is < 10 ng per mg antibody as determined by the LAL method.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	The antibody is purified using human lymphotoxin betaR affinity chromatography.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab10493** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

ELISA

WB

应用说明	ELISA: Use at a concentration of 0.5 - 1 µg/ml to detect a limit of ~0.16 ng/well of recombinant human lymphotoxin betaR. WB: Use at a concentration of 0.1 - 0.2 µg/ml to detect human lymphotoxin betaR at 5 ng/lane
------	---

and 25 ng/lane under non reducing and reducing conditions, respectively. Predicted molecular weight: 46 kDa.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

---

## 靶标

功能	Receptor for the heterotrimeric lymphotoxin containing LTA and LTB, and for TNFS14/LIGHT. Promotes apoptosis via TRAF3 and TRAF5. May play a role in the development of lymphoid organs.
序列相似性	Contains 4 TNFR-Cys repeats.
细胞定位	Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

---

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors