

## Product datasheet

# Anti-ZNF652 antibody ab105516

### 1 图像

#### 概述

---

产品名称	Anti-ZNF652抗体
描述	兔多克隆抗体to ZNF652
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Goat, Horse, Cow, Cat, Dog, Saccharomyces cerevisiae



免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 97-146 (RENSDDTEEE EEEVSYKREQ IVEVNLNNQ TLNVSKGEKG VSSQSKETPV) of Human ZNF652 (NP_055712). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
-----	--

阳性对照	HeLa cell lysate
------	------------------

#### 性能

---

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

---

Our [Abpromise guarantee](#) covers the use of **ab105516** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 0.5 µg/ml. Predicted molecular weight: 70 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

功能	Functions as a transcriptional repressor.
组织特异性	Widely expressed with higher expression in breast, prostate, vulva and pancreas.
序列相似性	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 9 C2H2-type zinc fingers.
细胞定位	Nucleus.

## 图片



Anti-ZNF652 antibody (ab105516) at 0.5 µg/ml + HeLa cell lysate at 10 µg

**Predicted band size:** 70 kDa

Gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors