

## Product datasheet

# Anti-DDDDK tag antibody [29E4.G7] - (Equivalent to FLAG antibodies from Sigma) (FITC) ab106221

### 概述

<b>产品名称</b>	Anti-DDDDK tag抗体[29E4.G7] - (Equivalent to FLAG抗体from Sigma) (FITC)
<b>描述</b>	小鼠单克隆抗体[29E4.G7] to DDDDK tag - (Equivalent to FLAG抗体from Sigma) (FITC)
<b>宿主</b>	Mouse
<b>偶联物</b>	FITC. Ex: 493nm, Em: 528nm
<b>特异性</b>	ab106221 recognizes the DYKDDDDK epitope tag (Asp-Tyr-Lys-Asp-Asp-Asp-Asp-Lys) fused to either the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.
<b>免疫原</b>	Synthetic peptide corresponding to DDDDK tag conjugated to Keyhole Limpet Haemocyanin (KLH).

### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.88% Sodium chloride, 1% BSA
<b>克隆</b>	单克隆
<b>克隆编号</b>	29E4.G7
<b>同种型</b>	IgG

### 靶标

**相关性** This is a useful tool for the localisation and characterisation of DDDDK tagged proteins.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors