

Product datasheet

Anti-IRGC antibody ab106483

3 图像

概述

产品名称	Anti-IRGC抗体
描述	兔多克隆抗体to IRGC
宿主	Rabbit
特异性	Two isoforms of IRGC are known to exist; ab106483 will recognize both isoforms.
经测试应用	适用于: ICC/IF, WB, ICC
种属反应性	与反应: Mouse, Rat, Human
免疫原	12 amino acid peptide near the carboxy terminus of Human IRGC (NP_062558)
阳性对照	Mouse brain tissue lysate A20 cells

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
存储溶液	Preservative: 0.02% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab106483** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
ICC/IF		Use a concentration of 20 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 50 kDa.
ICC		Use a concentration of 2.5 µg/ml.

## 靶标

### 序列相似性


Belongs to the interferon-inducible GTPase family.

## 图片



Immunofluorescence of IRGC in A20 cells with ab106483 at 20 ug/mL.

Immunocytochemistry/ Immunofluorescence - Anti-IRGC antibody (ab106483)



**Lane 1** : Anti-IRGC antibody (ab106483) at 1  $\mu$ g/ml

**Lane 2** : Anti-IRGC antibody (ab106483) at 2  $\mu$ g/ml

**All lanes** : Mouse brain tissue lysate

Lysates/proteins at 15  $\mu$ g per lane.

**Predicted band size:** 50 kDa

Western blot - IRGC antibody (ab106483)



Immunocytochemistry of IRGC in A20 cells  
with ab106483 at 2.5 µg/ml.

Immunocytochemistry - IRGC antibody (ab106483)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors