

## Product datasheet

# Anti-Endothelin 2 antibody ab106601

### 1 图像

#### 概述

产品名称	Anti-Endothelin 2抗体
描述	兔多克隆抗体to Endothelin 2
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MVSVPTTWCS VALALLVALH EGKGQAAATL EQPASSSHAQ GTHLRLRRCS) of Human Endothelin 2 (NP_001947). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	HepG2 cell lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab106601** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
----	------	----

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 20 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

功能	Endothelins are endothelium-derived vasoconstrictor peptides.
组织特异性	Expressed in lung, but not in placental stem villi vessels or cultured placental villi smooth muscle cells.
序列相似性	Belongs to the endothelin/sarafotoxin family.
细胞定位	Secreted.

## 图片



Anti-Endothelin 2 antibody (ab106601) at 1 µg/ml + HepG2 cell lysate at 10 µg

**Predicted band size:** 20 kDa

Gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors