

Product datasheet

Anti-P protein antibody ab106990

1 图像

概述

产品名称	Anti-P蛋白抗体
描述	兔多克隆抗体to P蛋白
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to N terminal amino acids 118-146 of Human P protein (NP_000266.2).
阳性对照	293 cell line lysates

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	This antibody is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab106990** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 93 kDa.

靶标

功能	Could be involved in the transport of tyrosine, the precursor to melanin synthesis, within the melanocyte. Regulates the pH of melanosome and the melanosome maturation. One of the components of the mammalian pigmentary system. Seems to regulate the post-translational processing of tyrosinase, which catalyzes the limiting reaction in melanin synthesis. May serve as a key control point at which ethnic skin color variation is determined. Major determinant of brown and/or blue eye color.
疾病相关	Defects in OCA2 are the cause of albinism oculocutaneous type 2 (OCA2) [MIM:203200]. An autosomal recessive disorder in which the biosynthesis of melanin pigment is reduced in skin, hair, and eyes. Although affected infants may appear at birth to have complete absence of melanin pigment, most patients acquire small amounts of pigment with age. Visual anomalies include decreased acuity and nystagmus. The phenotype is highly variable. The hair of affected individuals may turn darker with age, and pigmented nevi or freckles may be seen. African and African American individuals may have yellow hair and blue-gray or hazel irides. One phenotypic variant, 'brown OCA,' has been described in African and African American populations and is characterized by light brown hair and skin color and gray to tan irides.
序列相似性	Belongs to the CitM (TC 2.A.11) transporter family.
细胞定位	Melanosome membrane.

图片



Anti-P protein antibody (ab106990) at 1/100 dilution + 293 cell line lysates at 35 µg

Predicted band size: 93 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors