

Product datasheet

Anti-KCTD15 antibody ab107049

2 图像

概述

产品名称	Anti-KCTD15抗体
描述	兔多克隆抗体to KCTD15
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to C terminal amino acids 262-291 of Human KCTD15 (NP_076981.2, NP_001123466.1, NP_001123467.1).
阳性对照	K562 cell line lysates and mouse NIH-3T3 cell line lysates

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	<a href="#">ab10749</a> is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab107049** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 32 kDa.

靶标

## 序列相似性


Contains 1 BTB (POZ) domain.

## 图片



Anti-KCTD15 antibody (ab107049) at 1/100 dilution + K562 cell line lysates at 35  $\mu$ g

**Predicted band size:** 32 kDa



Anti-KCTD15 antibody (ab107049) at 1/100 dilution + Mouse NIH-3T3 cell line lysates at 35  $\mu$ g

**Predicted band size:** 32 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors