

Product datasheet

Anti-HAPLN4 antibody ab108156

1 图像

概述

产品名称	Anti-HAPLN4抗体
描述	兔多克隆抗体to HAPLN4
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Horse, Guinea pig, Cow, Cat, Dog
免疫原	Synthetic peptide corresponding to a region within internal sequence amino acids 251-300 (NYGYRHNAEE RYDAFCFTSN LPGRVFFLKP LRPVVPFSGAA RACAARGAAV) of Human HAPLN4 (NP_075378.1). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	HeLa cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab108156** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 43 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能	Binds to hyaluronic acid and may be involved in formation of the extracellular matrix.
组织特异性	Expressed predominantly in brain.
序列相似性	Belongs to the HAPLN family. Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 2 Link domains.
细胞定位	Secreted > extracellular space > extracellular matrix.

图片



Anti-HAPLN4 antibody (ab108156) at 1 µg/ml
+ HeLa cell lysate at 10 µg

Predicted band size: 43 kDa

Gel concentration 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors