

Product datasheet

Anti-STK32C antibody ab110453

2 图像

概述

产品名称	Anti-STK32C抗体
描述	兔多克隆抗体to STK32C
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide derived from internal sequence of Human STK32C.
阳性对照	Human breast carcinoma tissue; Extracts from Jurkat cells

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: PBS, 0.88% Sodium chloride, 50% Glycerol
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab110453** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 55 kDa.
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

靶标

序列相似性

Belongs to the protein kinase superfamily, Ser/Thr protein kinase family.
Contains 1 protein kinase domain.

图片




All lanes : Anti-STK32C antibody (ab110453)
at 1/500 dilution

Lane 1 : Extracts from Jurkat cells

Lane 2 : Extracts from Jurkat cells with
immunising peptide at 10 µg

Lysates/proteins at 30 µg per lane.

Predicted band size : 55 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - STK32C antibody (ab110453)

ab110453, at 1/50 dilution, staining STK32C
in paraffin-embedded Human breast
carcinoma tissue by Immunohistochemistry, in
the presence (right panel) or absence (left
panel) of immunising peptide.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors