

Product datasheet

Anti-B4GALT5 antibody ab110720

1 图像

概述

产品名称	Anti-B4GALT5抗体
描述	兔多克隆抗体to B4GALT5
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse 预测可用于: Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Pig, Zebrafish
免疫原	Synthetic peptide corresponding to a region within C-terminal amino acids 311-360 (GYSVSRPEGD TGKYKSIPHH HRGEVQFLGR YALLRKSKEK QGLDGLNNLN) of Mouse B4GALT5 (NP_062809). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Mouse Kidney Lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituents: PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab110720** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 45 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能	Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids.
组织特异性	Ubiquitously expressed.
通路	Protein modification; protein glycosylation.
序列相似性	Belongs to the glycosyltransferase 7 family.
细胞定位	Golgi apparatus > Golgi stack membrane. Trans cisternae of Golgi stack.

图片



Anti-B4GALT5 antibody (ab110720) at 1 $\mu\text{g}/\text{ml}$ + Mouse Kidney Lysate at 10 μg

Predicted band size: 45 kDa

12% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors