

Product datasheet

Anti-CD48 antibody [MRC OX-45] (Biotin) ab112160

1 References 1 图像

概述

产品名称	Anti-CD48抗体[MRC OX-45] (Biotin)
描述	小鼠单克隆抗体[MRC OX-45] to CD48 (Biotin)
宿主	Mouse
偶联物	Biotin
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Rat
免疫原	Rat T cell blasts (stimulated purified T helper cells with allogeneic irradiated rat spleen cells).
阳性对照	Rat Strain: Wister (Bone Marrow, Thymus)

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituents: 99% PBS, 0.5% BSA
纯度	Protein G purified
克隆	单克隆
克隆编号	MRC OX-45
骨髓瘤	NS0/1
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab112160** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		Use 0.2µg for 10 <sup>6</sup> cells. <a href="#">ab18434</a> -Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

## 靶标

功能	Ligand for CD2. Might facilitate interaction between activated lymphocytes. Probably involved in regulating T-cell activation.
序列相似性	Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Cell membrane.

## 图片



Cell Source: Thymus Percentage of cells stained above control: 94.9%. [ab112160](#) at a concentration of 0.2 µg.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors